# I931 Ionomètre de paillasse



Caractéristiques :

Ionomètre professionnel pour connecter toute électrode à ion spécifique Conçu pour une utilisation sur paillasse en laboratoire Ecran géant, Bras flexible porte électrode, alimentation sur secteur. 2 à 5 Points de Calibration : (0.001, 0.01, 0.1, 1, 10, 100, 1000, 10000ppm, mg/L, mol/L, mmol/L)

Lecture directe de la concentration en ion.

Différentes unités de mesure au choix : ppm, mg/L ou mol/L. La compensation automatique de température fournit une lecture précise sur

toute la gamme et une icône de stabilité affiche automatiquement l'état de la mesure du courant.

Fonction Hold pour geler la lecture sur l'écran et une mesure en mV est disponible pour afficher la performance de l'électrode

Diagnostic automatique de l'électrode avec affichage de la pente. Guide d'utilisation et d'aide en lecture sur l'écran géant de l'instrument Menu système entièrement paramètrable sur 8 critères avec choix du nombre de points de calibration, de la précision, de la stabilité de lecture, de l'arrêt automatique etc. Mémoire interne avec capacité de stockage de 500 données. Possibilité de transfert des données sur PC via un logiciel dédié.

#### Inclus:

Instrument, Bras flexible porte-électrode, Sonde de Temperature , cable USB et adaptateur prise secteur

Electrode (s) à ion spécifique et standards de calibration à choisir par l'utilisateur en fonction de ses besoins analytiques

Model	I931
Gamme Ion :	0.001~19999ppm, mg/L, mol/L (en fonction de la gamme de mesure de l'électrode à ion spécifique)
Précision Ion :	$\pm$ 0.5% pleine échelle (Monovalent), $\pm$ 1% pleine échelle (Divalent)
Gamme mV :	-1999.9~1999.9mV
Précision mV :	$\pm$ 0.2mV
Gamme Température :	0~105°C, 32~221°F
Précision Température :	±0.5°C, ±0.9°F
Points de Calibration :	2~5 points (0.001, 0.01, 0.1, 1, 10, 100, 1000, 10000ppm, mg/L, mol/L, mmol/L)
Temperature Compensation	0~100°C, 32~212°F, Manuelle ou Automatique
Fonction Hold :	Manuelle ou Automatique
Réglage de la stabilité	Basse ou Haute
Rappel calibration :	0 à 31 jours
Mémoire interne :	Stocke jusqu'à 500 données
Sortie :	Interface USB
Alimentation :	DC5V, avec adaptateur AC , 220VAC/50Hz
Dimmensions :	210(L)×188(W)×60(H)mm
Poids :	1.5kg



WATER KITS SUPPLY Tél:+33(0) 5 62 95 17 94 Email: contact@water-kits.fr Téléport 5 - 62590 Juillan - France

Web : www.water-kits.fr

### Introduction

Thank you for selecting the I931 benchtop ion meter. This manual provides a step-by-step guide to help you operate the meter, please carefully read the following instructions before use.

#### Unpacking the Meter

Before unpacking, ensure that the current work environment meet following conditions.

- Relative humidity is less than 80%.
- Ambient temperature is greater than 0°C and less than 60°C.
- No potential electromagnetic interference.

The following list describes the standard accessories of the meter. After the unpacking, please check all accessories are complete. If any are damaged or missing, please contact nearest distributor.

ACCESSORIES:

- Temperature Probe
- USB Cable
- DC5V Power Adapter

# Display

1931 benchtop ion meter is equipped with a clear and bright LCD display that used to show measured values, mode indicators and help message. The following table describes the meaning of each indicator.



#### INDEX:

Measure	Measurement mode indicator: Indicates meter is in the measurement mode	Stable	Stable indicator: Indicates the measured value has stabilized
Calibration	Calibration mode indicator: Indicates meter is in the calibration mode	Hold	Hold indicator: Indicates the displayed value has been frozen
Setup	Setup mode indicator: Indicates meter is in SETUP mode	Ð	Calibration Due Reminder: Prompts user to calibrate the meter regularly
Memory	Memory mode indicator: Indicates data is stored into memory	ATC	Automatic Temperature Compensation: Indicates meter is in the temperature compensation mode

# Keypad

The meter has a succinct membrane keypad, names and symbols describe the each function key controls.



#### Connectors

1931 benchtop ion meter provides 4 connectors for connecting the various types of sensors. Listed in the below table are the details of these connectors.



#### INDEX:

NO.	CONNECTOR	FUNCTION
1	ISE	For connecting the Ion Selective Electrode
2	ATC	For connecting the temperature probe
3	REF	For connecting the reference electrode
4	USB/ ம	For connecting the USB cable and DC5V power adapter

# Connecting the Electrode Holder

The meter comes with an easy-to-use holder for mounting a variety of Ion Selective Electrodes and temperature probe. If necessary, please follow the steps below to install the electrode holder.

1. The base of the electrode holder with an irregular round hole.



2. The electrode arm has a plastic rod.



3. Insert the plastic rod into the irregular round hole and swivel the electrode arm 90°. Electrode holder is now ready to swing into desired position.



#### ADJUSTMENT OF ELECTRODE ARM:

After installation, if the electrode arm automatically rises or falls, you need to adjust the screws until arm locate at any position.

1. Remove the plastic cover from the electrode arm.



2. Use the screwdriver to tighten the screw moderately.

3. Insert the plastic cover to previous position. Installation is completed.

#### Connecting the Electrode

1. Take out the Ion Selective Electrode from the packaging. Follow the steps below to place the electrode into left or right sides of the electrode arm.



 Insert the BNC connector into the corresponding connector socket. Rotate and push the connector clockwise until it locks. After connection is completed, DO NOT pull on the sensor cord. Always make sure that the connector is clean and dry.



### Prior to Use

- Remove the protective cap from the bottom of the Ion Selective Electrode.
- Soak the electrode in the diluted standard solution (e.g., 100ppm) for at least 10 minutes.



### Power On/Off

- Press MEAS key to turn on the meter, the display shows measured values, mode indicators and help messages.
- Press and hold the MEAS key for 3 seconds, the meter will turn off.
- If you do not press any key within the specified time period, the meter will turn off automatically.

# $(\mathbf{\hat{I}})$

To disable the auto-off function, please read the SETUP MENU section.

# Setup Menu

1931 benchtop ion meter contains an integrated setup menu that allows you to customize each displayed option to meet measurement requirements.

PARAMETER	DESCRIPTION	OPTIONS	DESCRIPTION	DEFAULT
		ppm	Parts per million	•
	Measurement, Unit:	mg/L	Milligrams per liter	
ПЦ ІЕ	Sets the default ion concentration and temperature	mol/L	Moles per liter	
	units.	°۲	Degrees Celsius	•
		°F	Degrees Fahrenheit	
		2	2 points	•
C 01	Calibration Points:	3	3 points	
LAL	Select the number of calibration points.	Ч	4 points	
		5	5 points	
100	Ion Valence:	1	Monovalent	•
1011	Select the ion valence of sensor.	2	Divalent	
528	Stable Criteria: Sets the stability criteria for measurement. When the "LO" option is enabled, measuring value	LO	Low	•
	will stabilize quickly, but reading is less accurate. When the "HI" option is enabled, measuring value will stabilize slowly, but guarantees high accuracy.	н	High	
801 년	Auto-Hold: When the auto-hold function is enabled, the meter	962	Enable	
,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,	will automatically sense a stable end-point reading and freeze it.	по	Disable	•
	Auto Dowon Off	10	10 minutes	
000	When the auto-off power is enabled, if you do not	20	20 minutes	
UFF	press any key within a specified time period, the meter will automatically turn off.	30	30 minutes	
		по	Disable	•
C 81 1	Calibration Due: When calibration due reminder is enabled, if you do	131	1 to 31 days	
	not recalibrate meter within a specified time period, the meter will automatically show $$ indicator.	OFF	Disable	•

#### 1931 Benchtop Ion Meter

98FE	Date and Time: Sets the date and time of the meter.			
<u> </u>	Clear stored data:	YES	Enable	
	Clear all stored data.	по	Disable	•
- 6 4	Reset: Reset function allows user to restore the meter back to factory default parameters. When this function is enabled, all calibration values and selected parameters will be lost or reset.	YE 5	Enable	
r 3C		по	Disable	•

SETTING THE DEFAULT PARAMETERS:

1. Press and hold the 🗎 key for 3 seconds, the meter enters setup menu, the display shows selectable parameter and page number.



2. Press 🛦 or 🔻 key to scroll through menu, select the parameter you want to set (Refer to Setup Menu section).

ION	Setup
ŀ	1019
	P-05
MESSAGE:	
Press ▲ or ▼ to select option or value     Press ENTER to confirm     Press MEAS to return to measurement	mode

3. Press ENTER key to confirm, the display shows an option in the submenu.

8



4. Press  $\blacktriangle$  or  $\blacktriangledown$  key to select the desired option.

5. Press ENTER key to confirm, the meter returns to measurement mode. Setting is completed.

#### EXIT THE SETUP MENU:

During the setup mode, if you want to exit setup menu, press MEAS key, the meter will return to measurement mode immediately.

### Setting the Date and Time

1931 benchtop ion meter has a real time clock that is used to time-stamp stored measured value and calibration data. Follow the steps below to set the date and time during the first use.

- 1. Press and hold the  $\, \boxminus \,$  key for 3 seconds to enter the setup menu.
- 2. Press  $\blacktriangle$  or  $\blacktriangledown$  key until the display shows "Date" option.



3. Press ENTER key to confirm, the meter shows current year.

ION	Setuo	
5	51.0	  Year
	YER-	
MESSAGE:		
Press ▲ or ▼ to select option or value     Press ENTER to confirm     Press MEAS to return to measurement mode		

- 4. Press  $\blacktriangle$  or  $\blacktriangledown$  key to set the year.
- 5. Press ENTER key to confirm, the meter shows current date and time ((Format: mm-dd, hh-mm)

ION	Setup	
	730-	Format: Month - Day
	12:55	Format: Hour - Minutes
MESSAGE:		-
Press ▲ or ▼ to select option or value     Press ENTER to confirm     Press MEAS to return to measurement mode		

7. Press ENTER key to confirm, the meter returns to measurement mode. Setting is completed.

### Selecting the Calibration Points

1931 benchtop ion meter supports ion concentration calibration up to 5 points with minimum of 2 points, available calibration points include the following options.

MEASUREMENT UNITS	CALIBRATION POINTS
ppm	0.001, 0.01, 0.1, 1, 10, 100, 1000, 10000
mg/L	0.001, 0.01, 0.1, 1, 10, 100, 1000, 10000
mol/L	0.001, 0.01, 0.1, 1, 10
mmol/L	0.001, 0.01, 0.1

If you need to modify the number of calibration points, please follow the steps below.

1. Press and hold the  $\square$  key for 3 seconds to enter the setup menu.

2. Press A key, the display shows "CAL/P-02" (Calibration Point) option.



3. Press ENTER key to confirm, the meter enters the setting mode.

4. Press  $\blacktriangle$  or  $\blacktriangledown$  key to select the number of calibration points.

ION	Setup
	2
	ERL
MESSAGE:	
<ul> <li>Press ▲ or ▼ to select option or value</li> <li>Press ENTER to confirm</li> </ul>	
· Press MEAS to return to measurement mode	

5. Press ENTER key to confirm, the meter returns to measurement mode. Setting is completed.

### Selecting the Ion Valence

The meter is capable of connecting a variety of ion selective electrodes. For the divalent ions, you need to set the ion valence before calibration or measurement.

- 1. Press and hold the  $\, \boxplus \,$  key for 3 seconds to enter the setup menu.
- 2. Press 🔺 key, the display shows "ION/P-03" (Ion Valence) option



- 3. Press ENTER key to confirm, the meter enters the setting mode.
- 4. Press  $\blacktriangle$  or  $\blacktriangledown$  key to select the ion valence.

2
חסו

5. Press ENTER key to confirm, the meter returns to measurement mode. Setting is completed.

#### Selecting the Concentration Unit

I931 benchtop ion meter is capable of using the mg/L, ppm or mol/L as measurement unit of concentration. The factory default is ppm. If you need to convert measurment unit, the meter must be recalibrated.

- 1. Press and hold the  $\hfill \hfill \hfill$
- 2. Press  $\blacktriangle$  key, the display shows "UNIT" option.



- 3. Press ENTER key to confirm, the meter enters the setting mode.
- 4. Press  $\blacktriangle$  or  $\blacktriangledown$  key to select the desired concentration unit (ppm, mg/L, mol/L).
- 5. Press ENTER key to confirm, the "CAL" indicator will flashing uninterruptedly indicating that the meter is waiting for calibrating.

ION	Setup	
	ERL	mgiL
MESSAGE:		
Press CAL to enter calibration mode		
Press MEAS to return to measurement mode		

6. Press CAL key to enter the calibration mode or MEAS key to cancel option.

### **Temperature Compensation**

In order to get accurate measuring results, you need to enable the manual or automatic temperature compensation before measurement or calibration.

AUTOMATIC TEMPERATURE COMPENSATION:

• Place the temperature probe into the circular hole of electrode arm.



• Insert the connector of temperature probe into corresponding connector socket (Marked "ATC").



• The "ATC" indicator will show on the display, the meter is now switched to automatic temperature compensation mode.



MANUAL TEMPERATURE COMPENSATION:

- 1. DO NOT connect the temperature probe to meter.
- 2. Press and hold the °C key for 3 seconds to enter temperature setting mode.
- 3. Press  $\blacktriangle$  or  $\blacktriangledown$  key to set the temperature value of sample.
- 4. Press ENTER key to confirm, the meter returns to measurement mode. Setting is completed.

# $(\mathbf{\hat{U}})$

In the temperature setting mode, press  $\blacktriangle$  or  $\checkmark$  key once, the setting value will increase or decrease by 0.1. Press and hold the  $\bigstar$  or  $\checkmark$  key, the setting value will increase or decrease by 1.

### Ion Concentration Calibration

To obtain accurate measurement results, we recommend that you perform ion calibration and measurement at same temperature. If you are not calibrate the meter or calibration is not successfully, the display will always show "0.000".

1. Press MODE key until the meter shows **ION** indicator.



2. Press CAL key, the meter shows 0.001ppm (or mg/L, mol/L).

		Calibration	
	0.		ppm
		ERL	1
MESSAGE:			
<ul> <li>Immerse the sensor</li> <li>Press ▲ or ▼ to se</li> <li>Press ENTER to co</li> </ul>	in the calibration solution lect option or value nfirm		
Press MEAS to retu	rn to measurement mode		

3. If necessary, press  $\blacktriangle$  or  $\checkmark$  key to select the desired calibration point (e.g., 0.01ppm).

ION	Calibration	
	0.0	ppm
	ERL	1
MESSAGE:		
· Immerse the senser is the collication solution		
<ul> <li>Press A or V to select option or value</li> <li>Press ENTER to confirm</li> </ul>		
Press MEAS to return to measurement mode		

- 4. Rinse the Ion Selective Electrode with distilled water, then rinse with a small amount of ion standard solution.
- 5. Dip the electrode into corresponding calibration solution. Stir the sensor gently to create a homogenous solution.



6. Press ENTER key to confirm, "Calibration" indicator begins flashing. Wait for the measured value to stabilize, the display shows "0.1/CAL2". The meter prompts you to continue with second point calibration.

- 7. Rinse the Ion Selective Electrode with distilled water. Dip the electrode into corresponding calibration solution. Stir the sensor gently.
- 8. Press ENTER key, "Calibration" indicator begins flashing. Wait for the measured value to stabilize, the display shows "1/CAL3". The meter prompts you to continue with third point calibration.



9. Repeat steps 7 to 8 above until the display shows "END", the meter returns to measurement mode automatically. Calibration is completed.



#### EXIT THE CALIBRATION:

During the calibration process, if you want to exit calibration, press MEAS key, the meter will return to measurement mode immediately.

### Ion Calibration Report

This program lets you check the slope of the Ion Selective Electrode.

- 1. Press MR key in the ion measurement mode, the meter shows "LOC/P-01".
- 2. Press  $\blacktriangle$  or  $\triangledown$  key until the display shows "ELE/P-02" (Electrode Diagnosis).



3. Press ENTER key to confirm, the meter shows the last calibration date (Format: mm-dd).

ION	Memory
	2-03
MESSAGE:	
Press ▲ or ▼ to select option or value     Press MEAS to return to measurement mode	

4. Press  $\mathbf{\nabla}$  key, the meter shows calibration point and its slope.



5. After the browsing, press MEAS key to exit the current mode.

### **Temperature Calibration**

During the measurement, when automatic temperature compensation is enabled, if the temperature reading displayed differs from that of an accurate thermometer, you need to calibrate the meter.

1. Press and hold the °C key for 3 seconds to enter temperature calibration mode, the display shows current temperature reading.



- 2. Press  $\blacktriangle$  or  $\blacktriangledown$  key to set the temperature value.
- 3. Press ENTER key to confirm. Calibration is completed.

#### Ion Concentration Measurement

- 1. Press MODE key until the meter shows **ION** indicator.
- 2. Rinse the Ion Selective Electrode thoroughly with distilled water. Dip the electrode into the sample solution.
- 3. If your sample is belong to low concentration liquids or some interfering ions are present in the solution, we suggest you that adding the lonic Strength Adjuster into the sample solution.
- 4. Stir the sensor gently. Wait for the reading to stabilize, record the measured value on the display.

#### mV Measurement

Press MODE key until display shows measurement unit "mV". Rinse the electrode thoroughly with distilled water. Dip the electrode into the sample solution. Wait for the measured value to stabilize, record the reading on the display.

#### Hold Function

1931 benchtop ion meter contains two data hold modes. When the Auto-Hold function is enabled, the meter will automatically sense a stable endpoint reading and freeze it, "HOLD" indicator appears on the display. If the Auto-Hold function is disabled, press  $\widehat{\bullet}$  key, the meter will immediately freeze currently displayed value. Press the key again to resume measuring.



### Storing and Recalling Data from Memory

The meter allows up to 500 data sets to be stored and recalled.

#### MEMORY INPUT:

During the measurement process, press MI key to input measured value into the memory, "Memory" indicator appears on the display.



MEMORY RECALL:

- 1. Press MR key in the measurement mode, the meter shows "LOC" (Data Log).
- 2. Press ENTER key to confirm, the meter shows page number of the stored data.

ION	Memory	
	10	Page number
MESSAGE:		
<ul> <li>Press ▲ or ▼ to select option or value</li> </ul>		
Press MEAS to return to measurement mode		

3. Press  $\mathbf{\nabla}$  key, the meter shows date and time of the stored data (Format: mm-dd).



4. Press  $\mathbf{\nabla}$  key again, the display shows the stored data.



5. After the browsing, press MEAS key to exit the current mode.

### Ion Selective Electrode Care and Maintenance

- Ensure that the electrode is thoroughly washed with distilled water after each use.
- DO NOT scratch the sensitive membrane on electrode during the measurement or calibration.
- If you do not use the electrode for long periods, please store the electrode in a dry, cool and well-ventilated area.

# Troubleshooting

LCD DISPLAY	CAUSE	CORRECTIVE ACTION
	Electrode dried out	Soak the Ion Selective Electrode in the diluted standard solution at least 15 minutes
	Measured value is out of range	Check the electrode whether clogged, dirty or broken
Err	Incorrect calibration solutions	Using the fresh calibration solutions for calibration
	Electrode is out of service life	Replace the Ion Selective Electrode

# Specifications

lon	Model	1931	
	Range	0.001~19999ppm, mg/L, mol/L (Depending on range of ISE)	
	Accuracy	±0.5% F.S (Monovalent), ±1% F.S (Divalent)	
	Calibration Points	2 to 5 points	
	Calibration Solutions	0.001, 0.01, 0.1, 1, 10, 100, 1000, 10000ppm, mol/L, mg/L	
mV	Range	-1999.9~1999.9mV	
	Accuracy	±0.2mV	
	Resolution	0.1mV	
Temperature	Range	0~105°C, 32~221°F	
	Accuracy	±0.5°C, ±0.9°F	
	Resolution	0.1°C	
	Calibration Points	1 point, Measured value $\pm 10^{\circ}$ C	
Others	Temperature Compensation	0~100°C, 32~212°F, Manual or Automatic	
	Hold Function	Manual or Automatic	
	Stability Conditions	Low or High	
	Calibration Due	O to 31 days	
	Power Off	Manual or Automatic (10, 20, 30 minutes)	
	Memory	Stores up to 500 data sets	
	Output	USB Communication Interface	
	Power Requirements	DC5V, using AC adapters, 220VAC/50Hz	
	Dimensions	210(L)×188(W)×60(H)mm	
	Weight	1.5kg	